

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD24

Numéro de catalogue: **CL647-65121**



Informations de base

Numéro de catalogue:

CL647-65121

Taille:

100tests , 5 µl/test

Hôte:

Rat

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

BC075622

Identification du gène (NCBI):

12484

Nom complet:

CD24a antigen

Méthode de purification:

Purification par affinité

CloneNo.:

M1/69

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and 5.00 ul CoraLite® Plus 647-conjugated Anti-Mouse CD24 (CL647-65121, Clone: M1/69). Cells were not fixed.



1X10⁶ C57BL/6 mouse splenocytes were surface stained with FITC-Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and CoraLite® Plus 647-conjugated rat IgG2b isotype control antibody. Cells were not fixed.