

À des fins de recherche uniquement

Anticorps Monoclonal anti-TER-119

Numéro de catalogue: CL647-65149



Informations de base

Numéro de catalogue:

CL647-65149

Taille:

100tests , 5 µl/test

Hôte:

Rat

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

104231

Nom complet:

Lymphocyte antigen 76

Méthode de purification:

Purification par affinité

CloneNo.:

TER-119

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

TER-119 is a mouse erythroid lineage-specific monoclonal antibody that reacts with erythroid cells at differentiation stages from early proerythroblast to mature erythrocyte, but not with cells showing typical erythroid blast-forming unit (BFU-E) and erythroid colony-forming unit (CFU-E) activities (PMID: 1975515; 10848813). TER-119 recognizes a 52-kDa molecule on erythrocyte membranes (PMID: 10848813). TER-119 antigen is a molecule associated with cell-surface glycoporphin A but not with glycoporphin A itself (PMID: 10848813).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ C57BL/6 mouse bone marrow cells were surface stained with FITC-Anti-Mouse CD45 (FITC-65087, Clone: 30-F11) and 5.00 ul CoraLite® Plus 647-conjugated Anti-Mouse TER-119 (CL647-65149, Clone: TER-119). Cells were not fixed.



1X10⁶ C57BL/6 mouse bone marrow cells were surface stained with FITC-Anti-Mouse CD45 (FITC-65087, Clone: 30-F11) and CoraLite® Plus 647-conjugated Rat IgG2b isotype control antibody. Cells were not fixed.