

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD62P

Numéro de catalogue: **CL647-65164**



Informations de base

Numéro de catalogue:

CL647-65164

Taille:

100tests , 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC028067

Identification du gène (NCBI):

6403

Nom complet:

selectin P (granule membrane protein
140kDa, antigen CD62)

MW calculé

830 aa, 91 kDa

Méthode de purification:

Purifié par chromatographie d'affinité
par protéine A

CloneNo.:

AK4

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain, primates non humains

Informations générales

P-selectin, also known as GMP140 or CD62P, is a transmembrane glycoprotein that mediates the interaction of activated endothelial cells or platelets with leukocytes. It is an adhesion molecule involved in the pathogenesis of inflammation, thrombosis, and oncogenesis. P-selectin is stored in the alpha-granules of platelets and Weibel-Palade bodies of endothelial cells. Upon cell activation by agonists, P-selectin is transported rapidly to the cell surface.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Thrombin-activated human peripheral blood platelets were surface co-stained with FITC Anti-Human CD41a and 5 ul CoraLite® Plus 647 Anti-Human CD62P (CL647-65164, Clone: AK4) or unstained. Cells were not fixed.