

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD58

Numéro de catalogue: **CL647-65193**



Informations de base

Numéro de catalogue:

CL647-65193

Taille:

100tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC005930

Identification du gène (NCBI):

965

Nom complet:

CD58 molecule

MW calculé

28 kDa

Méthode de purification:

Purification par protéine G

CloneNo.:

TS2/9

Excitation/Emission maxima
wavelengths:

654 nm / 674 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 5 ul CoraLite® Plus 647 Anti-Human CD58 (CL647-65193, Clone:TS2/9) or unstained. Cells were not fixed. Monocytes were gated.