

À des fins de recherche uniquement

Anticorps Monoclonal anti-Integrin beta-7



Numéro de catalogue: **CL647-65213**

Informations de base

Numéro de catalogue:
CL647-65213

Taille:
100tests , 5 µl/test

Hôte:
Rat

Isotype:
IgG2a, kappa

Numéro d'acquisition GenBank:
BC011184

Identification du gène (NCBI):
16421

Nom complet:
integrin beta 7

Méthode de purification:
Purification par protéine A

CloneNo.:
FIB504

Excitation/Emission maxima
wavelengths:
654 nm / 674 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
Humain, souris

Informations générales

Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS avec azoture de sodium à 0,1 % et BSA à 0,5 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface co-stained with Anti-Human CD3 and 5 ul CoraLite® Plus 647 Anti-Human/Mouse Integrin beta-7 (CL647-65213, Clone: FIB504), or Rat IgG2a Isotype Control. Cells were not fixed. Lymphocytes were gated.