

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD71

Numéro de catalogue: **CL647-65235**



## Informations de base

Numéro de catalogue:

CL647-65235

Taille:

100tests , 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC001188

Identification du gène (NCBI):

7037

Nom complet:

transferrin receptor (p90, CD71)

MW calculé

85 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

OKT9

Excitation/Emission maxima  
wavelengths:

654 nm / 674 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



$1 \times 10^6$  PHA-activated human PBMCs were surface stained with 5  $\mu$ l CoralLite® Plus 647 Anti-Human CD71 (CL647-65235, Clone: OKT9), or APC Mouse IgG1 Isotype Control. Cells were not fixed.