

À des fins de recherche uniquement

# Anticorps Monoclonal anti-IgG2a Isotype Control



Numéro de catalogue: **CL647-65244**

## Informations de base

Numéro de catalogue:

CL647-65244

Taille:

100ug, 0.2 mg/ml

Hôte:

Mouse

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

MOPC-173

Excitation/Emission maxima  
wavelengths:

654 nm / 674 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

## Informations générales

## Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug CoraLite® Plus 647 Anti-Mouse IgG2a Isotype Control (CL647-65244, Clone: MOPC-173) (red), CoraLite® Plus 647-Conjugated Anti-Human CD8 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.