

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD49b

Numéro de catalogue: **CL647-65245**



Informations de base

Numéro de catalogue:
CL647-65245

Taille:
100ug , 500 µg/ml

Hôte:
Armenian Hamster

Isotype:
IgG

Numéro d'acquisition GenBank:
BC065139

Identification du gène (NCBI):
16398

Nom complet:
integrin alpha 2

Méthode de purification:
Purification par affinité

CloneNo.:
HMA2

Excitation/Emission maxima
wavelengths:
654 nm / 674 nm

Applications

Applications testées:
FC

Spécificité de l'espèce:
souris

Informations générales


Stockage

Stockage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Tampon de stockage:
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse NK1.1 (CD161) (FITC-65138, Clone: PK136) and 0.5 ug CoraLite® Plus 647 Anti-Mouse CD49b (CL647-65245, Clone: HMA 2) or 0.5 ug CoraLite® Plus 647-conjugated Armenian Hamster IgG. Cells were not fixed.