

À des fins de recherche uniquement

Anticorps Monoclonal anti-TRAF6

Numéro de catalogue: **CL647-66498**

Phare



Informations de base

Numéro de catalogue:

CL647-66498

Taille:

100ul, Concentration: 1000 µg/ml by Nanodrop;

Hôte:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24696

Numéro d'acquisition GenBank:

BC031052

Identification du gène (NCBI):

7189

Nom complet:

TNF receptor-associated factor 6

MW calculé

522 aa, 60 kDa

MW observés:

60 kDa

Méthode de purification:

Purification par protéine G

CloneNo.:

1D1E1

Excitation/Emission maxima wavelengths:

654 nm / 674 nm

Applications

Applications testées:

FC (Intra)

Spécificité de l'espèce:

Humain

Informations générales

TRAF6 is a member of the TNF receptor associated factor (TRAF) protein family. TRAF6 is recognized as a signal transducer, which activates the NF-κB pathway in response to pro-inflammatory cytokines. TRAF6 has E3 ligase activity and is responsible for inducing Lys-63 (K63)-linked poly-ubiquitination chains, functioning together with E2 Ubc13/Uev1A complex to mediate IKK activation. TRAF6 plays a crucial role in the regulation of the innate and adaptive immune responses.

Stockage

Stockage:

Stocker à -20 °C. Éviter toute exposition à la lumière. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec glycérol à 50 %, Proclin300 à 0,05 % et BSA à 0,5 %, pH 7,3.

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ Jurkat cells were intracellularly stained with 0.2 ug Coralite® Plus 647 Anti-Human TRAF6 (CL647-66498, Clone:1D1E1) (red), or 0.2 ug APC-65128 ; APC Mouse IgG2b Isotype Control (APC-65128, Clone: MPC-11) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).