

À des fins de recherche uniquement

Anticorps Monoclonal anti-PLVAP

Numéro de catalogue: **CL750-65214**



Informations de base

Numéro de catalogue:

CL750-65214

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC013517

Identification du gène (NCBI):

84094

Nom complet:

plasmalemma vesicle associated protein

Méthode de purification:

Purification par protéine A

CloneNo.:

MECA-32

Excitation/Emission maxima wavelengths:

755 nm / 780 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ bEnd.3 cells were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse PLVAP (CL750-65214, Clone:MECA-32) (red), or 0.5 ug CoraLite® Plus 750 Rat IgG2a Isotype Control (2A3) (CL750-65209, Clone: 2A3) (blue). Cells were not fixed.



1X10⁶ bEnd.3 cells were surface stained with 0.5 ug CoraLite® Plus 750 Anti-Mouse PLVAP (CL750-65214, Clone:MECA-32), or 0.5 ug CoraLite® Plus 750 Rat IgG2a Isotype Control (2A3) (CL750-65209, Clone: 2A3). Cells were not fixed.