

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD95

Numéro de catalogue: FITC-65071



Informations de base

Numéro de catalogue:

FITC-65071

Taille:

100 tests, 10 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC012479

Identification du gène (NCBI):

355

Nom complet:

Fas (TNF receptor superfamily, member 6)

MW calculé

35-38 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

DX2

Excitation/Emission maxima wavelengths:

494 nm / 520 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Fas (CD95/APO-1) is a transmembrane glycoprotein belonging to the tumor necrosis factor (TNF) receptor superfamily. It can mediate apoptosis by ligation with an agonistic anti-Fas antibody or Fas ligand. Stimulation of Fas results in the aggregation of its intracellular death domains, leading to the formation of the death-inducing signaling complex (DISC). FAS-mediated apoptosis may have a role in the induction of peripheral tolerance, in the antigen-stimulated suicide of mature T-cells, or both.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 %.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human peripheral blood lymphocytes were surface stained with 10 ul FITC Anti-Human CD95 (FITC-65071, Clone:DX2) (green) or FITC-mouse IgG1 isotype control (black). Cells were not fixed .