

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD4

Numéro de catalogue: FITC-65134

1 Publications



## Informations de base

Numéro de catalogue:

FITC-65134

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

BC025782

Identification du gène (NCBI):

920

Nom complet:

CD4 molecule

MW calculé

55 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

OKT4

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

## Applications

Applications testées:

FC

Demandes citées:

FC

Spécificité de l'espèce:

Humain

## Informations générales

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

## Publications notables

Autrice	Pubmed ID	Journal	Application
Guopin Pan	35131482	Pharmacol Res	FC

## Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA, pH 7.3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> human PBMCs were surface co-stained with Anti-Human CD3 and 5 ul FITC Plus Anti-Human CD4 (FITC-65134, Clone:OKT4), or 5 ul Mouse IgG2b Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul FITC Plus Anti-Human CD4 (FITC-65134, Clone:OKT4), or 5 ul Mouse IgG2b Isotype Control. Cells were not fixed. Lymphocytes were gated.