

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD4

Numéro de catalogue: FITC-65147



Informations de base

Numéro de catalogue:

FITC-65147

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC025782

Identification du gène (NCBI):

920

Nom complet:

CD4 molecule

MW calculé

55 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

SK3

Excitation/Emission maxima
wavelengths:

494 nm / 520 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

Tampon à base de phosphate avec azoture de sodium à 0,09 % et gélatine à 0,1 %, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD3 (APC-65151, clone UCHT1) and 5.00 ul/test FITC Anti-Human CD4 (FITC-65147, Clone:SK3). Sample were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD3 (APC-65151, clone UCHT1) and FITC mouse IgG1 isotype control antibody. Sample were not fixed.