

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD1a

Numéro de catalogue: FITC-65172



Informations de base

Numéro de catalogue:

FITC-65172

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC031645

Identification du gène (NCBI):

909

Nom complet:

CD1a molecule

MW calculé

327 aa, 37 kDa

Méthode de purification:

L'anticorps purifié est conjugué à l'isothiocyanate de fluorescéine (FITC) dans des conditions optimales. Le conjugué est purifié par chromatographie d'exclusion de taille.

CloneNo.:

HI149

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ MOLT-4 cells were surface stained with 5 ul FITC Plus Anti-Human CD1a (FITC-65172, Clone:HI149) or Mouse IgG1 Isotype Control. Cells were not fixed.