

À des fins de recherche uniquement

Anticorps Monoclonal anti-STRO-1

Numéro de catalogue: FITC-65184

1 Publications



Informations de base

Numéro de catalogue:

FITC-65184

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgM, lambda

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purifié par précipitation et chromatographie

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

Applications

Applications testées:

FC

Demandes citées:

FC

Spécificité de l'espèce:

Humain

Espèces citées:

Humain

Informations générales

Publications notables

Autrice	Pubmed ID	Journal	Application
Tao Jiang	37441553	J Tissue Eng	FC

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human whole blood cells were surface stained with 5 ul FITC Plus Anti-Human STRO-1 (FITC-65184, Clone: STRO-1) or unstained. Cells were not fixed. CD45dim, CD34+ cells were gated.