

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG2a Isotype Control



Numéro de catalogue: FITC-65208

Informations de base

Numéro de catalogue:

FITC-65208

Taille:

100ug, 0.5 mg/mL

Hôte:

Mouse

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

C1.18.4

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.2 ug FITC Anti-Mouse IgG2a Isotype Control (FITC-65208, Clone:C1.18.4) (red) or 0.2 ug FITC Anti-Mouse CD45.2 (FITC-65072) (green), or unstained (blue). Cells were not fixed.