

À des fins de recherche uniquement

# Anticorps Monoclonal anti-PLVAP

Numéro de catalogue: FITC-65214



## Informations de base

Numéro de catalogue:

FITC-65214

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC013517

Identification du gène (NCBI):

84094

Nom complet:

plasmalemma vesicle associated protein

Méthode de purification:

Purification par protéine A

CloneNo.:

MECA-32

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> bEnd.3 cells were surface stained with 1 ug FITC Plus Anti-Mouse PLVAP (FITC-65214, Clone: MECA-32) (red) or 1 ug Control Antibody. Cells were not fixed.