

À des fins de recherche uniquement

# Anticorps Monoclonal anti-IgM Isotype Control



Numéro de catalogue: FITC-65243

## Informations de base

Numéro de catalogue:

FITC-65243

Taille:

100ug , 200 µg/ml

Hôte:

Mouse

Isotype:

IgM, lambda

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

N/A

CloneNo.:

MOPC-104E

Excitation/Emission maxima  
wavelengths:

495 nm / 524 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

## Informations générales

## Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> human whole blood cells were surface stained with 0.2 ug FITC Plus Anti-Mouse IgM Isotype Control (red), 5 ul FITC Plus Anti-Human STRO-1 (green), or unstained (blue). Cells were not fixed.