

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD223

Numéro de catalogue: PE-65098



Informations de base

Numéro de catalogue:

PE-65098

Taille:

100ug, 0.2 mg/ml

Hôte:

Rat

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC120591

Identification du gène (NCBI):

16768

Nom complet:

Lymphocyte-activation gene 3

Méthode de purification:

Purification par affinité

CloneNo.:

C9B7W

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Lymphocyte activation gene 3 protein (LAG-3), also known as CD223, is a type I transmembrane protein belonging to the immunoglobulin superfamily (PMID: 1692078). It is a CD4-related molecule that binds MHC class II with high affinity (PMID: 1692078; 9159144). LAG-3 is expressed on activated T and NK cells. Besides, natural CD4(+)CD25(+) Tregs express LAG-3 upon activation (PMID: 15485628). LAG-3 negatively regulates T-cell function and homeostasis, and contributes to the suppressor activity of Tregs (PMID: 14644131; 15485628).

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.1% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ anti-CD3/CD28 stimulated (3 days) mouse splenocytes were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD3 and 0.2 ug PE Anti-Mouse CD223 (PE-65098, Clone: C9B7W) or 0.2 ug PE Mouse IgG1 Isotype Control. Cells were not fixed.