

À des fins de recherche uniquement

Anticorps Monoclonal anti-LPAM-1

Numéro de catalogue: PE-65101



Informations de base

Numéro de catalogue:

PE-65101

Taille:

100ug, 0.2 mg/ml

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC011184

Identification du gène (NCBI):

16421

Nom complet:

integrin beta 7

Méthode de purification:

Purification par affinité

CloneNo.:

DATK32

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse bone marrow cells were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD45R (B220) (RA3-6B2) and 0.2 ug PE Anti-Mouse LPAM-1 (PE-65101, Clone:DATK32) or 0.2 ug Isotype Control. Cells were not fixed.