

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD24

Numéro de catalogue: PE-65121



Informations de base

Numéro de catalogue:

PE-65121

Taille:

100ug, 0.2 mg/ml

Hôte:

Rat

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

BC075622

Identification du gène (NCBI):

12484

Nom complet:

CD24a antigen

Méthode de purification:

Purification par affinité

CloneNo.:

M1/69

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

Tampon à base de phosphate avec azoture de sodium à 0,09 % et gélatine à 0,1 %, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.06 ug PE Anti-Mouse CD24 (PE-65121, Clone: M1/69) (red) or 0.06 ug isotype control antibody (black). Cells were not fixed.