

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD62P

Numéro de catalogue: PE-65164



## Informations de base

Numéro de catalogue:

PE-65164

Taille:

100 tests, 5 µl/test

Hôte:

Mouse

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC028067

Identification du gène (NCBI):

6403

Nom complet:

selectin P (granule membrane protein  
140kDa, antigen CD62)

MW calculé

830 aa, 91 kDa

Méthode de purification:

L'anticorps purifié est conjugué à la R-phycoérythrine (PE) dans des conditions optimales. Le conjugué est purifié par chromatographie d'exclusion de taille.

CloneNo.:

AK4

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain, primates non humains

## Informations générales

P-selectin, also known as GMP140 or CD62P, is a transmembrane glycoprotein that mediates the interaction of activated endothelial cells or platelets with leukocytes. It is an adhesion molecule involved in the pathogenesis of inflammation, thrombosis, and oncogenesis. P-selectin is stored in the alpha-granules of platelets and Weibel-Palade bodies of endothelial cells. Upon cell activation by agonists, P-selectin is transported rapidly to the cell surface.

## Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> thrombin-activated human peripheral blood platelets were surface co-stained with Anti-Human CD41a labeled with FlexAble CoraLite® Plus750 Antibody Labeling Kit for Mouse IgG1 (KFA024) and PE Anti-Human CD62P (PE-65164, Clone: AK4) or Mouse IgG1 Isotype Control. Cells were not fixed.