

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG2a Isotype Control



Numéro de catalogue: PE-65209

Informations de base

Numéro de catalogue:

PE-65209

Taille:

100ug, 0.2 mg/mL

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

2A3

Excitation/Emission maxima
wavelengths:

496 nm, 565 nm / 578 nm

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Stockage

Stockage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Tampon de stockage:

Tampon à base de phosphate avec azoture de sodium à 0,09 % et gélatine à 0,1 %, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.5 ug PE Rat IgG2a Isotype Control (PE-65209, Clone: 2A3) (blue) or 0.5 ug PE Anti-Mouse CD4 (PE-65141, Clone: RM4-5) (red). Cells were not fixed.