


RAS RECEPTOR PATHWAY

www.ptglab.com



Related Products

| Legend | |
|--------|--------------------|
| | Association |
| | Directly Activates |
| | Inhibits |
| | Phosphorylates |

| Antibody Name | Catalog Number | Type | Applications |
|---------------|----------------|-------------|----------------------------|
| AKT/PKB | 10176-2-AP | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| AP1, JUN, P59 | 10024-2-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| ATF2-Specific | 14834-3-AP | Rabbit Poly | ELISA, IHC, WB |
| ATF6 | 24169-3-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| CD3 | 17617-3-AP | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| CDC42 | 10155-3-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| CDK1 | 10762-1-AP | Rabbit Poly | ELISA, WB |
| CREB1 | 12208-1-AP | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| ERK1/2 | 16443-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| GRB2 | 10254-2-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| GSR/GR | 18257-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| HRAS-Specific | 18295-1-AP | Rabbit Poly | ELISA, FC, IF, IHC, WB |
| IKKβ/IKK1 | 15649-1-AP | Rabbit Poly | ELISA, IF, WB |

| Antibody Name | Catalog Number | Type | Applications |
|---------------|----------------|-------------|----------------------------|
| JNK | 10023-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| KRAS | 12063-1-AP | Rabbit Poly | ELISA, FC, IF, IHC, IP, WB |
| KRAS | 60309-1-ig | Mouse Mono | ELISA, FC, IHC, WB |
| MYC | 10828-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| NFκB1 | 14220-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| NRAS | 10724-1-AP | Rabbit Poly | ELISA, FC, IF, IHC, WB |
| p38 MAPK | 14064-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| PKC1 | 10754-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PCNA/Cyclin | 10205-2-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| PIK3CA/PI3K | 20583-1-AP | Rabbit Poly | ELISA, FC, IHC, WB |
| PLD1 | 18355-1-AP | Rabbit Poly | ELISA, IF, IP, WB |
| PLK3/CNK | 10977-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| PTEN | 60300-1-ig | Mouse Mono | ELISA, IF, IHC, WB |

| Antibody Name | Catalog Number | Type | Applications |
|---------------|----------------|-------------|------------------------|
| RAB5A | 11947-3-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| RAC1 | 24072-3-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| RAF1 | 51140-1-AP | Rabbit Poly | ELISA, IHC, WB |
| RALA/RAL | 13629-1-AP | Rabbit Poly | ELISA, IHC, WB |
| RAP1B | 10840-1-AP | Rabbit Poly | ELISA, WB |
| RASA1/RASGAP | 12955-1-AP | Rabbit Poly | ELISA, IHC, IP, WB |
| RASGRF1 | 12958-1-AP | Rabbit Poly | ELISA, IP, WB |
| RHOA | 10749-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| RRAF1 | 11647-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| SYK | 12206-1-AP | Rabbit Poly | ELISA, IHC, IP, WB |
| TEK/TIE2 | 19157-1-AP | Rabbit Poly | ELISA, IF, IHC, WB |
| VAV1 | 16364-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |
| ZAP70 | 15592-1-AP | Rabbit Poly | ELISA, IF, IHC, IP, WB |

This number shows the amount of times our antibody has been cited in a publication.