


p38-MAPK SIGNALING

www.ptglab.com



Related Products

Legend	
	Association
	Directly Activates
	Inhibits
	Binding
	Phosphorylates

Antibody Name	Catalog Number	Type	Applications
ASK1	14385-1-AP	Rabbit Poly	ELISA, WB
ATF2-Specific	14834-3-AP	Rabbit Poly	ELISA, IHC, WB
Caspase 3	19677-3-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CDC25A	55031-3-AP	Rabbit Poly	ELISA, IHC, IP, WB
CDC25B	10644-3-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
CDC25C	16485-3-AP	Rabbit Poly	ELISA, IP, WB
CDC42	10155-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
EGFR (C-terminal)	51071-2-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB
EGFR (N-terminal)	22542-1-AP	Rabbit Poly	ELISA, FC, IHC, WB
IL-1 alpha	16765-1-AP	Rabbit Poly	ELISA, IF, WB
IL-1 beta	16806-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
IRAK 1	10478-2-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
JUN	66313-1-Ig	Mouse Mono	ELISA, WB

Antibody Name	Catalog Number	Type	Applications
MEK2	55106-1-AP	Rabbit Poly	ELISA, IP, WB
MEK3	21072-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
MNK1	10136-1-AP	Rabbit Poly	ELISA, IF, WB
MSK1	15977-1-AP	Rabbit Poly	ELISA, IF, WB
NFKB1	23576-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
NFKB2	51085-1-AP	Rabbit Poly	ELISA, WB
p38 MAPK	66234-1-Ig	Mouse Mono	ELISA, IHC, WB
pan Ras	60309-1-Ig	Mouse Mono	ELISA, FC, IHC, WB
PKC alpha	66421-1-Ig	Mouse Mono	ELISA, IF, IHC, WB
PKC beta	132919-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
PKC gamma	14364-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, WB
Rac1	24072-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
RAC1/2/3	16117-1-AP	Rabbit Poly	ELISA, IF, IHC, IP, WB

Antibody Name	Catalog Number	Type	Applications
RAC2	10735-1-AP	Rabbit Poly	ELISA, IHC, IP, WB
SAPK4	10217-1-AP	Rabbit Poly	ELISA, IHC, WB
TAB1	14819-1-AP	Rabbit Poly	ELISA, IP, WB
TAB2	14410-1-AP	Rabbit Poly	ELISA, FC, IF, IHC, IP, WB
TAK1	26250-2-AP	Rabbit Poly	ELISA, IHC, WB
TAOK1	26250-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
TAOK2	21188-1-AP	Rabbit Poly	ELISA, IP, WB
TAOK3	10158-2-AP	Rabbit Poly	ELISA, IHC, IP, WB
TGF-beta 1	18978-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
TGF-beta 2-Specific	19999-1-AP	Rabbit Poly	ELISA, IF, IHC, WB
TGF-beta 3	18942-1-AP	Rabbit Poly	ELISA, IHC, WB
TRAF6	12809-1-AP	Rabbit Poly	ELISA, IF, WB

This number shows the amount of times our antibody has been cited in a publication.