




Speedy One-Step ELISA Kits

Proteintech's solution for rapid analyte detection in a broad range of sample types

- Results in 90 minutes or less
- Reduced assay complexity with single-wash step protocols
- Unparalleled precision and sensitivity
- FcZero™ Capture Antibodies
- Unmatched Cost savings



Scan to Learn More



Speedy ELISA kits are made with antibody pairs containing capture antibodies with Proteintech's FcZero backbone. These antibodies have been engineered to silence the Fc-receptor binding site, thus eliminating non-specific binding with Fc receptors sometimes present in serum and plasma samples.

Visit ptglab.com for more information




Speedy ELISA Kit Offering

| Target | Catalog | Sensitivity |
|--------------------|---------|-------------|
| Human APOE | SE50006 | 0.03 ng/mL |
| Human HMGB1 | SE50009 | 0.05 ng/mL |
| Human AFP | SE50010 | 7.32 ng/mL |
| Human TNF-alpha | SE50011 | 2.81 ng/mL |
| Human OX40L/TNFSF4 | SE50012 | 0.03 ng/mL |
| Human IL-6 | SE50013 | 1.06 ng/mL |
| Mouse IL-2 | SE51002 | 1.71 ng/mL |
| Rat IL-2 | SE52002 | 2.67 ng/mL |

Speedy Kits enable analyte detection in 90 minutes or less. That's nearly 3x faster than a traditional ELISA! Simply add all prepared samples, standards, and antibodies in a single incubation step, wash, develop, and read.

Rapid, Single-Wash Step Protocol



Visit ptglab.com for more information