



chromotek® DYKDDDDK Fab-Trap®

A premium resin for high purity immunoprecipitation of FLAG®-tagged proteins.

- Lowest background in class
- Less contaminating antibody chains
- Recombinantly expressed Fab-fragment provides reproducible results
- Compatible with harsh buffer conditions
- Ready-to-use



▲ DYKDDDDK Fab-Trap® consisting of an anti-FLAG® Fab-fragment (blue) conjugated to Agarose beads (grey).



▲ Immunoprecipitation of FLAG® fusion protein with DYKDDDDK Fab-Trap®, elution with 3xDYKDDDDK-peptide. I: Input, FT: Flow-Through, E: Elution.

FLAG® is a registered trademark of Sigma-Aldrich Co LLC.

chromotek[®] DYKDDDDK reagents for Immunoprecipitation & Western blot

- DYKDDDDK Fab-Trap[®] Agarose
- 3xDYKDDDDK-peptide
- DYKDDDDK tag Polyclonal antibody
- DYKDDDDK Immunoprecipitation Starter Pack



**Scan to order your
free sample**

For more information visit, ptglab.com