

For Research Use Only

BCRP/ABCG2 Polyclonal antibody

Catalog Number: 10051-1-AP **14 Publications**



Basic Information

Catalog Number:

10051-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0005

GenBank Accession Number:

BC021281

GeneID (NCBI):

9429

UNIPROT ID:

Q9UNQ0

Full Name:

ATP-binding cassette, sub-family G (WHITE), member 2

Calculated MW:

70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IHC

Species Specificity:

human

Cited Species:

human, mouse, canine, chicken

Notable Publications

Author	Pubmed ID	Journal	Application
Jingsheng Yuan	31078735	Cancer Lett	WB
Manxin Xie	29567542	Biomed Pharmacother	WB
Ping Xie	28569775	Cell Death Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data