

For Research Use Only

CCNI Polyclonal antibody

Catalog Number: 10279-1-AP



Basic Information

Catalog Number: 10279-1-AP	GenBank Accession Number: BC000420	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 547 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 10983	Recommended Dilutions: IF 1:50-1:500
Source: Rabbit	Full Name: cyclin I	
Isotype: IgG	Calculated MW: 43 kDa	
Immunogen Catalog Number: AG0300		

Applications

Tested Applications: IF, ELISA	Positive Controls: IF : HepG2 cells,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data