

For Research Use Only

Goat-IgG Polyclonal antibody

Catalog Number: 10285-1-AP

50 Publications



Basic Information

Catalog Number:

10285-1-AP

GenBank Accession Number:

GeneID (NCBI):

Purification Method:

Antigen affinity purification

Size:

150ul, Concentration: 1400 ug/ml by Nanodrop and 400 ug/ml by Bradford method using BSA as the standard;

Full Name:

Calculated MW:

155 kDa

Source:

Rabbit

Isotype:

IgG

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF, RIP

Species Specificity:

goat

Cited Species:

mouse

Notable Publications

Author	Pubmed ID	Journal	Application
Jian Wen	32996777	Cancer Biother Radiopharm	WB
Lan Zhang	34586738	Clin Transl Med	RIP
Hucheng Chen	30272275	Mol Med Rep	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1%BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

