

For Research Use Only

Calmodulin 1/2/3 Polyclonal antibody



Catalog Number: 10303-1-AP

Basic Information

Catalog Number: 10303-1-AP	GenBank Accession Number: BC000454	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 100 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 801	
Source: Rabbit	Full Name: calmodulin 1 (phosphorylase kinase, delta)	
Isotype: IgG	Calculated MW: 17 kDa	
Immunogen Catalog Number: AG0490		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human lung cancer using 10303-1-AP (CALM1 antibody) at dilution of 1:100 (under 25x lens).