

For Research Use Only

ZNF256 Polyclonal antibody, PBS Only

Catalog Number: 10360-1-PBS



Basic Information

Catalog Number:

10360-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0390

GenBank Accession Number:

BC001438

GeneID (NCBI):

10172

UNIPROT ID:

Q9Y2P7

Full Name:

zinc finger protein 256

Calculated MW:

55 kDa

Observed MW:

72 kDa, 55 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

ZNF256 as a transcriptional repressor that plays a role in cell proliferation. ZNF256 localizes in the nucleus and contains the KRAB domain that mediates interaction with TRIM28 and is required for transcriptional repressor activity. ZNF256 exists as two isoforms and the calculated molecular weight of isoforms is 55 and 75 kDa.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 10360-1-AP (ZNF256 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 10360-1-PBS in a different storage buffer formulation.