

MMP15 Polyclonal antibody

Catalog Number: 10372-1-AP

Basic Information

Catalog Number:

10372-1-AP

Size:

150ul , Concentration: 247 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0548

GenBank Accession Number:

BC006412

GeneID (NCBI):

4324

UNIPROT ID:

P51511

Full Name:

matrix metalloproteinase 15 (membrane-inserted)

Calculated MW:

72 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Background Information

MMP15(Matrix metalloproteinase-15) is a endopeptidase that degrades various components of the extracellular matrix.MMP15-specific antibodies detect a 72-kD protein in lung cell surface membranes, suggesting that MMP15 is membrane-bound. They demonstrated by Northern blotting that MMP15 is widely expressed as a 3.6-kb transcript, with highest levels in liver, placenta, testis, colon, and intestine(PMID:7649159).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20° C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

