

For Research Use Only

# Calmodulin Polyclonal antibody

Catalog Number: 10521-1-AP **2 Publications**



## Basic Information

<b>Catalog Number:</b> 10521-1-AP	<b>GenBank Accession Number:</b> BC006464	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 253 µg/ml by Bradford method using BSA as the standard;	<b>GeneID (NCBI):</b> 805	
<b>Source:</b> Rabbit	<b>Full Name:</b> calmodulin 2 (phosphorylase kinase, delta)	
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 17 kDa	
<b>Immunogen Catalog Number:</b> AG0820		

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB

**Species Specificity:**  
human

**Cited Species:**  
mouse

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Heng-Xue Wang	34484665	Toxicol Res (Camb)	WB
Jennie Sjöhamn	27643892	Protein Sci	

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data