

For Research Use Only

BL21(DE3)-lysate Polyclonal antibody

Catalog Number: 10558-1-Ig **1 Publications**



Basic Information

Catalog Number:

10558-1-Ig

GenBank Accession Number:

GeneID (NCBI):

Size:

150ul, Concentration: 2000 ug/ml by Nanodrop and 700 ug/ml by Bradford method using BSA as the standard;

Full Name:

Calculated MW:

1-300 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:1000-1:4000

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0048

Applications

Tested Applications:

WB, ELISA

Positive Controls:

WB : Total bacteria lysates of BL21(DE3),

Cited Applications:

WB

Species Specificity:

bacteria

Cited Species:

zebrafish

Notable Publications

Author	Pubmed ID	Journal	Application
Sakura Tanaka	34978328	Endocrinology	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

