

For Research Use Only

BL21(DE3)-lysate Polyclonal antibody



Catalog Number: 10558-1-Ig

1 Publications

Basic Information

Catalog Number: 10558-1-Ig	GenBank Accession Number: GeneID (NCBI):	Purification Method: Protein A purification
Size: 150ul , Concentration: 2000 µg/ml by Nanodrop and 700 µg/ml by Bradford method using BSA as the standard;	Full Name: Calculated MW: 1-300 kDa	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit		
Isotype: IgG		
Immunogen Catalog Number: AG0048		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Total bacteria lysates of BL21(DE3),
Cited Applications: WB	
Species Specificity: bacteria	
Cited Species: zebrafish	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Sakura Tanaka	34978328	Endocrinology	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

