

For Research Use Only

# His-Bacteria-Protein Polyclonal antibody



Catalog Number: 10560-1-Ig **1 Publications**

## Basic Information

Catalog Number:

10560-1-Ig

Size:

150ul , Concentration: 773 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0050

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

30 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:2000-1:16000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : 30 kDa His-tag bacteria protein of the E.Coli(DE3),

## Notable Publications

Author	Pubmed ID	Journal	Application
Ziling Fang	29547122	Elife	

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

