

For Research Use Only

# BL21(plysS) -lysate Polyclonal antibody

Catalog Number:10562-1-Ig



## Basic Information

Catalog Number:

10562-1-Ig

Size:

150ul , Concentration: 2000 ug/ml by Nanodrop and 700 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0045

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Calculated MW:

1-200 kDa

Purification Method:

Protein A purification

Recommended Dilutions:

WB 1:1000-1:4000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

bacteria

Positive Controls:

WB : Total bacteria lysates of BL21(plysS),

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

