

For Research Use Only

# ECSIT Polyclonal antibody

Catalog Number: 10577-1-AP



---

## Basic Information

|  |   |  |
|--|---|--|
| <b>Catalog Number:</b><br>10577-1-AP   | <b>GenBank Accession Number:</b><br>BC008279    | <b>Purification Method:</b><br>Antigen affinity purification |
| <b>Size:</b><br>150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard; | <b>GeneID (NCBI):</b><br>51295                  |  |
| <b>Source:</b><br>Rabbit   | <b>Full Name:</b><br>ECSIT homolog (Drosophila) |  |
| <b>Isotype:</b><br>IgG   | <b>Calculated MW:</b><br>49 kDa                 |  |
| <b>Immunogen Catalog Number:</b><br>AG0850   |   |  |

---

## Applications

**Tested Applications:**  
IHC, ELISA

**Species Specificity:**  
human

---

## Background Information

---

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


---

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon cancer using 10577-1-AP (ECSIT antibody) at dilution of 1:50 (under 10x lens).