

For Research Use Only

PTGES Polyclonal antibody

Catalog Number: 10596-1-AP



Basic Information

Catalog Number:

10596-1-AP

Size:

150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0847

GenBank Accession Number:

BC008280

GeneID (NCBI):

9536

UNIPROT ID:

O14684

Full Name:

prostaglandin E synthase

Calculated MW:

17 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

Human prostaglandin E synthase (PTGES) is a member of a recently recognized protein superfamily consisting of membrane associated proteins involved in eicosanoid and glutathione metabolism (the MAPEG family) (PMID:10377395). It is also named as MGS1L1, MPGES1, PGES, PIG12 and belongs to the MAPEG family. This protein is expressed in membranes as an approximately 16-kD protein (PMID:10091672).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

