

For Research Use Only

Endothelin 3 Polyclonal antibody

Catalog Number: 10674-1-AP **1 Publications**



Basic Information

Catalog Number: 10674-1-AP	GenBank Accession Number: BC008876	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 133 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1908	
Source: Rabbit	Full Name: endothelin 3	
Isotype: IgG	Calculated MW: 26 kDa	
Immunogen Catalog Number: AG1071		

Applications

Tested Applications:
IHC, ELISA

Cited Applications:
IF

Species Specificity:
human, mouse, rat

Cited Species:
mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Juliano T Freitas	34288510	Pigment Cell Melanoma Res	IF

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 10674-1-AP (Endothelin 3 antibody) at dilution of 1:100 (under 10x lens).