

For Research Use Only

# Endothelin 3 Polyclonal antibody

Catalog Number: 10674-1-AP

1 Publications



## Basic Information

Catalog Number:

10674-1-AP

Size:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1071

GenBank Accession Number:

BC008876

GeneID (NCBI):

1908

UNIPROT ID:

P14138

Full Name:

endothelin 3

Calculated MW:

26 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

IHC, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse

## Notable Publications

Author	Pubmed ID	Journal	Application
Juliano T Freitas	34288510	Pigment Cell Melanoma Res	IF

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human breast cancer using 10674-1-AP (Endothelin 3 antibody) at dilution of 1:100 (under 10x lens).