

For Research Use Only

SUFU Polyclonal antibody

Catalog Number: 10836-1-AP

2 Publications



Basic Information

Catalog Number:

10836-1-AP

Size:

150ul , Concentration: 250 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1275

GenBank Accession Number:

BC013291

GeneID (NCBI):

51684

UNIPROT ID:

Q9UMX1

Full Name:

suppressor of fused homolog (Drosophila)

Calculated MW:

54 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Yasunori Sato	30036521	Am J Pathol	IF
Chia-Hao Huang	27377153	Lab Chip	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

