

For Research Use Only

CRK Polyclonal antibody

Catalog Number: 10973-1-AP



Basic Information

Catalog Number: 10973-1-AP	GenBank Accession Number: BC009837	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 1398	
Source: Rabbit	Full Name: v-crk sarcoma virus CT10 oncogene homolog (avian)	
Isotype: IgG	Calculated MW: 40 kDa	
Immunogen Catalog Number: AG1420		

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human prostate cancer using 10973-1-AP (CRK antibody) at dilution of 1:50 (under 10x lens).