

For Research Use Only

CENTD3 Polyclonal antibody

Catalog Number: 11024-3-AP



Basic Information

Catalog Number: 11024-3-AP	GenBank Accession Number: BC009968	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 167 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 64411	
Source: Rabbit	Full Name: ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3	
Isotype: IgG	Calculated MW: 8 kDa	
Immunogen Catalog Number: AG1479		

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

This gene encodes a phosphoinositide binding protein containing ARF-GAP, RHO-GAP, RAS-associating, and pleckstrin homology domains. The ARF-GAP and RHO-GAP domains cooperate in mediating rearrangements in the cell cytoskeleton and cell shape. It is a specific PtdIns(3,4,5)P3/PtdIns(3,4)P2-stimulated Arf6-GAP protein.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data