

For Research Use Only

LUZP4 Polyclonal antibody

Catalog Number: 11072-1-AP



Basic Information

Catalog Number: 11072-1-AP	GenBank Accession Number: BC128134	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 500 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 51213	Recommended Dilutions: WB 1:200-1:500
Source: Rabbit	Full Name: leucine zipper protein 4	
Isotype: IgG	Calculated MW: 40 kDa	
Immunogen Catalog Number: AG1539	Observed MW: 36 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : MDA-MB-231 cells, SKOV-3 cells, U-87 MG cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 11072-1-AP (LUZP4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.