

For Research Use Only

NOTCH2 Polyclonal antibody

Catalog Number: 11197-1-AP



Basic Information

Catalog Number: 11197-1-AP	GenBank Accession Number: BC019835	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 388677	
Source: Rabbit	Full Name: Notch homolog 2 (Drosophila) N-terminal like	
Isotype: IgG	Calculated MW: 265 kDa	
Immunogen Catalog Number: AG1679		

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Notch2 is another homologue gene of Norch in Drosophila. It is a 300kd transmembrane protein, with one 110kd processed fragment, on western blotting. Notch2, like other Norch genes, functions in development of various tissues. This antibody can recognize NOTCH2 and NOTCH2NL.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data