

For Research Use Only

# NOTCH2/NOTCH2NL Polyclonal antibody



Catalog Number: 11197-1-AP

1 Publications

## Basic Information

### Catalog Number:

11197-1-AP

### Size:

150ul , Concentration: 147 µg/ml by Bradford method using BSA as the standard;

### Source:

Rabbit

### Isotype:

IgG

### Immunogen Catalog Number:

AG1679

### GenBank Accession Number:

BC019835

### GeneID (NCBI):

388677

### UNIPROT ID:

PODPK4

### Full Name:

Notch homolog 2 (Drosophila) N-terminal like

### Calculated MW:

265 kDa

### Purification Method:

Antigen affinity purification

## Applications

### Tested Applications:

ELISA

### Cited Applications:

WB, IHC

### Species Specificity:

human

### Cited Species:

human

## Background Information

Notch2 is another homologue gene of Norch in Drosophila. It is a 300kd transmembrane protein, with one 110kd processed fragment, on western blotting. Notch2, like other Norch genes, functions in development of various tissues. This antibody can recognize NOTCH2 and NOTCH2NL.

## Notable Publications

Author	Pubmed ID	Journal	Application
Hansen Lin	38462037	Cancer Lett	WB,IHC

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

