

For Research Use Only

RBMX Polyclonal antibody

Catalog Number: 11205-1-AP



Basic Information

Catalog Number:

11205-1-AP

Size:

150ul , Concentration: 233 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1659

GenBank Accession Number:

BC007435

GeneID (NCBI):

27316

UNIPROT ID:

P38159

Full Name:

RNA binding motif protein, X-linked

Calculated MW:

43 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

RBMX, also named as HNRPG, RBMXP1 and Glycoprotein p43, is a rRNA-binding protein which may be involved in pre-mRNA splicing. Investigations into the RBM gene family are uncovering networks of protein interactions which regulate RNA processing, and which might operate downstream of signal transduction pathways.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

