

For Research Use Only

PRMT2 Polyclonal antibody

Catalog Number: 11269-1-AP



Basic Information

Catalog Number: 11269-1-AP	GenBank Accession Number: BC000727	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 140 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 3275	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: protein arginine methyltransferase 2	
Isotype: IgG	Calculated MW: 49 kDa	
Immunogen Catalog Number: AG1806		

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: human, mouse, rat	

Background Information

PRMT2 (also known as HRMT1L1) is a member of the arginine methyltransferase family which catalyzes arginine methylation and regulates diverse cellular processes, including transcription, translation, DNA repair, and RNA processing. Recently PRMT2 has been identified as a pro-tumorigenic factor in malignant gliomas.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data